

# Domtest project – outline and status

Instrumentation workshop LBNL

July 23 – 26 2003

TestDAQ Section

Martin Kestel for the

Penn State University DomTest group:

David Atlee, Jacob Carpenter, Phil Roth,

Doug Rutledge, Doug Cowen

# Outline


---

- DOM testing sequence:  
Interaction of DAQ-Control and analysis
- How domtests work
- More about tools

# A fixed sequence is repeated for each temperature step.


---

Ramp temperature down, measure (I, II, III)  
and repeat...



# In the freezer data recording is always on.


---



- I: Set target Temperature, wait for stabilisation
- II: Analyse data from section I, apply result. Determine nominal HV, apply.
- III: All other measurements with nominal settings.

# In part I the discriminator threshold is determined.


---



- Set pulser to various amplitudes, including the target SPE value.
- Determine the discriminator threshold and the gain for each DOM.

# In part I also the pedestal pattern is examined.

---



- Collect pedestal patterns at all temperatures
- Study gain per channel  
(see also Kael's talk)

At the end of part I DAQ Control asks  
the analysis package for the results.

---

- DAQ Control starts an analysis of the files recorded in part I.
- Results: discriminator threshold and pedestal patterns are immediately applied.
- Next: need nominal HV to continue → Part II

In part II the nominal HV is determined.

---

- DAQ Control starts a gain-versus-HV run.
- DAQ Control asks the analysis package to determine the nominal HV for the current temperature.
- How this is done → Following slide.



# DAQ Control asks for nominal HV settings for each DOM.

---

- DC specifies all filenames and output directory for the analysis package (xml).
- DC starts the analysis, passing xml as argument.
- Analysis writes (xml) results to output directory.
- DC reads in results and adjusts descriptor xmls for each DOM.
- DC asks all DOMs to re-adjust.

# In part III DAQ Control carries out a bunch more tests.

---

- Linearity tests using LEDs and a reference DOM
  - PMT time jitter
  - PMT quantum efficiency
  - Time synchronisation resolution (RapCal)
  - Many more....
- 
- Data get recorded but are analysed only later.

# Outline

---

- DOM testing sequence:  
Interaction of DAQ-Control and analysis
- How domtests work
- More about tools


# A domest checks a series of files for a specific property.

---

- Goal: use 10 weeks' data to judge DOM quality.
- Specify (xml):
  - Inputfiles
  - Suitable data subsets
  - Analysis class
  - Judgement class
  - Output directory
- Example:
  - /data/domtest/23July2003/TestDaq\_231.zip  
/data/domtest/23July2003/TestDaq\_232.zip
  - HV stable  
Temperature stable —▶ next slide  
Quality good
  - GainVersusHV —▶ later slide
  - GainVersusHVJudgement —▶ later slide
  - /data/domtest/23Jul2003/gainVsHV/results

# Identify time intervals where Temperature and HV were stable.






---



# Time intervals are specified e.g. as: 'HV was stable'.

Monitor Analysis Class Diagram

- Fixed limits:  
Value was in range [Low; Hi].
- Derived limits:  
Value was in a peak of time  
distribution of that quantity.



# The selected time intervals are stored in an xml file.

---

```
1 <?xml version="1.0" encoding="UTF-8"?>
2 <xs:schema xmlns:xs="http://www.w3.org/2001/XMLSchema">
3   <xs:include schemaLocation="MonitorRecordEntry.xsd"/>
4
5   <xs:element name="IntervalGroup">
6     <xs:complexType>
7       <xs:sequence>
8         <xs:element name="TimeIntervals" type="IntervalEntry" maxOccurs="unbounded"/>
9         <xs:element name="DACValues" type="DACEntry" maxOccurs="unbounded"/>
10      </xs:sequence>
11      <xs:attribute name="DomID" type="xs:string" use="required"/>
12    </xs:complexType>
13  </xs:element>
14
15  <xs:complexType name="IntervalEntry">
16    <xs:sequence>
17      <xs:element name="MonitorRecordEntry" type="MonitorEntry" maxOccurs="unbounded"/>
18    </xs:sequence>
19    <xs:attribute name="BeginTime" type="FortyEightBits" use="required"/>
20    <xs:attribute name="EndTime" type="FortyEightBits" use="required"/>
21  </xs:complexType>
22
23  <xs:complexType name="DACEntry">
24    <xs:attribute name="Name" type="xs:string" use="required"/>
25    <xs:attribute name="Value" type="xs:double" use="required"/>
26  </xs:complexType>
27
28 </xs:schema>
```

# Analysis is only carried out for records within valid time intervals.

---

- Book a set of result histograms for all independent time intervals.

Gain versus HV: find number of temperature intervals.

- Analyse all records from within time intervals.

Find pulse height and fill pulse-height-spectrum for the right temperature interval.

- Once the loop is finished, call 'finalise' method from the analysis.


Fit  $1/x$  and Gaussian to pulse-height-spectrum. Find peak position. This is a measure of the PMT gain @ HV (@ Temperature).

- Store the results.







# Analysis deals with waveforms extracted from domHitRecords.

Waveform Analysis Class Diagram




# Examples how the waveform analysis works – (1) baseline determination\*



(\*) Taken from the domtest documentation

# Examples how the waveform analysis works – (2) pulse height\*

Pulse Height Algorithm



(\*) Taken from the domtest documentation

# Analysis is only carried out for records within valid time intervals (2).

---

- Book a set of result histograms for all independent time intervals.

Gain versus HV: find number of temperature intervals.

- Analyse all records from within time intervals.

Find pulse height and fill pulse-height-spectrum for the right temperature interval.

- Once the loop is finished, call 'finalise' method from the analysis.


Fit  $1/x$  and Gaussian to pulse-height-spectrum. Find peak position. This is a measure of the PMT gain @ HV (@ Temperature).

- Store the results.

# DOM judgement for the gain versus HV measurement\*.

---

- Read in results from the analysis output.xml file.



- Judge e.g.  $\chi^2$  of distributions with straight lines and cut on their values: pass, fail, incomplete.

(\*) Not a terribly good example.

# Outline

---

- DOM testing sequence:  
Interaction of DAQ-Control and analysis
- How domtests work
- More about tools

# Other tools.


---

- Display package.
- Oscilloscope and Multi Channel Analyser with monitor value strip chart tool.
- Monitor strip chart as plugin for DAQ Control.
- DOM data simulator: test robustness of analysis with artificial glitches.

- Data reader/writer with data version awareness:

DomDescriptor (DomConfig), runHeader, domhitRecord, domMonitorRecord, extMonitorRecord, rapCalRecord.

# Display package offers help building up graphical displays.





# The display package provides a.o. a pre-factured main window.

---

## Constructors:

- TopLevelWindow()
- TopLevelWindow(String title, int length, int width, components, menus)

## Methods:

- addComponent(JComponent comp)
- addMenu(JMenu menu)
- displayWindow()

All other packages will use this main window and the logging interface.



# The display package provides also a pre-factured popup window.

---

## Constructors:

- PopupWindow()
- PopupWindow(String title, int length, int width, components)


## Methods:

- addComp(JComponent comp)
- setUpPanel()




The popup window is used when only a simple window is needed.

# The oscilloscope has four main parts.



# Simulator will be used to fake faulty data.

---



# The display package was used to create several tools (1).

---

Interface for  
waveform analysis.


The screenshot shows a software window titled "Waveform Analysis Test". It contains several panels and controls:

- Get waveforms**: A button at the top.
- Waveform Panel**: Contains a list box labeled "Waveforms" with indices 0 through 5.
- Waveform Calculation: 2**: A section with various parameters:
  - Peaks**: A text box containing the value "3".
  - Max Voltage Algorithm**: A dropdown menu set to "Three Level".
  - Baseline Algorithm**: A dropdown menu set to "Iterative".
  - Maximum Voltage**: A text box containing "126.74491610453889".
  - Baseline**: A text box containing "-0.2025592619073796".
  - Maximum Voltage Error**: A text box containing "5.196152422706632".
  - Baseline Error**: A text box containing "0.09662819047168711".
  - Pedestal RMS**: A text box containing "0.9516773076005155".
  - Pedestal RMS Error**: A text box containing "3.507253731803746E-5".
- Analyze Waveform**: A button below the calculation section.
- Peak Panel**: Contains a list box labeled "Peaks" with indices 0, 1, and 2.
- Peak Calculations**: A section with:
  - Peak Number**: A text box.
  - Max Voltage Algorithm**: A dropdown menu set to "Default".
- Log**: A section with a "Log" button and a "Log Level: OFF" label.
- Exit**: A button in the bottom right corner.

# The display package was used to create several tools (2).

---




Interface for  
data simulator.



# The display package was used to create several tools (3).

---

Interface to adjust DOM properties.



# Summary, outlook.

---

- Procedure setup in the freezer is converging
- A set of auxiliary tools (osci, mca, simulator) under way
- Versioned data formats under construction
- Release available soon
- Bottlenecks: reader/writer with versioning  
no hardware to play with (adds a lot of motivation)
- Need more input about: monitoring  
data basing procedures.